



Shanghai Joy Crown Industry Co., Ltd.

Add: Rm.1105 No.55-56 Lane 21 Pu San Road, Pu Dong New Area,

Shanghai 200125, P.R. China

Tel: +86-21-58751019 Fax: +86-21-50907607 Mobile/Whatsapp: 008613818337572

Email: sales@joycrown.com <https://www.medicalsupplyment.com>

Cohesive Fixation Bandage Specifications

Catalogue No. /Reference No. /Products Code	JC04026
Material	Fabric Ingredients: Approx. 40% Cotton and 60% Nylon, Coated with Latex Free Synthetic Adhesive
Weight per Sqm	75g±3g/m ²
Color	Brown
Size	Width: 4cm, 6cm, 8cm, 10cm, 12cm; Length (Stretched): 4m or 20m;
Elasticity	170%~190%
Latex or Latex Free	Latex Free
With/out Fastener	Without Clip, Velcro or other Fastener
Package	1 Roll/Shrink film
Sterile or Not	No Sterile
Rules & Classification	Rule 1, Class I





Shanghai Joy Crown Industry Co., Ltd.

Add: Rm.1105 No.55-56 Lane 21 Pu San Road, Pu Dong New Area,

Shanghai 200125, P.R. China

Tel: +86-21-58751019 Fax: +86-21-50907607 Mobile/Whatsapp: 008613818337572

Email: sales@joycrown.com <https://www.medicalsupplyment.com>

Cohesive Fixation Bandage Specifications

Reference No/ Product Code	Description	Size	Unit Package	Outside/Carton Package
JC040261	Cohesive Fixation Bandage, Different Sizes, 1 roll/Shrink film, non sterile.	2.5cm x 4m	1 Roll/Shrink film, 10 Rolls/Box	160 Rolls
JC040262		4cm x 4m	1 Roll/Shrink film, 5 Rolls/Box	160 Rolls
JC040263		6cm x 4m		160 Rolls
JC040264		8cm x 4m		80 Rolls
JC040265		10cm x 4m		80 Rolls
JC040266		12cm x 4m		80 Rolls
JC040267		4cm x 20m	1 Roll/Shrink film, 1 Rolls/Box	80 Rolls
JC040268		6cm x 20m		80 Rolls
JC040269		8cm x 20m		40 Rolls
JC0402610		10cm x 20m		40 Rolls
JC0402611		12cm x 20m		40 Rolls