



Shanghai Joy Crown Industry Co., Ltd.

Add: Rm.1105 No.55-56 Lane 21 Pu San Road, Pu Dong New Area,

Shanghai 200125, P.R. China

Tel: +86-21-58751019 Fax: +86-21-50907607 Mobile/Whatsapp: 008613818337572

Email: sales@joycrown.com <https://www.medicalsupplyment.com>

Tubular Net Bandage Specifications

Catalogue No. /Reference No. /Products Code	JC04015
Material	45% polyamide fiber; 55% latex
Weight per Sqm	≥ 260 gsm
Color	White
Size	0#, 0.5#, 1.0#, 2#, 3#, 4#, 5#, 5.5#, 5.8#, 6#, 7#, 8#, 9#, 10#.
Package	1 Roll/Polybag/Dispenser Box
Sterile or Not	No Sterile
Rules & Classification	Rule 1, Class I





Shanghai Joy Crown Industry Co., Ltd.

Add: Rm.1105 No.55-56 Lane 21 Pu San Road, Pu Dong New Area,

Shanghai 200125, P.R. China

Tel: +86-21-58751019 Fax: +86-21-50907607 Mobile/Whatsapp: 008613818337572

Email: sales@joycrown.com <https://www.medicalsupplyment.com>

Tubular Net Bandage Specifications

Reference No/ Product Code	Specification	Width	Stretched Length	Unit Package/Box	Outside/Carton Package
JC040151	0#	0.6cm	25m	1Roll/Polybag /Dispenser Box	180 box
JC040152	0.5#	0.9cm	25m		180 box
JC040153	1#	1.7cm	25m		120 box
JC040154	2#	2.0cm	25m		120 box
JC040155	3#	2.3cm	25m		84 box
JC040156	4#	2.7cm	25m		84 box
JC040157	5#	3.0cm	25m		84 box
JC040158	5.5#	3.5cm	25m		56 box
JC040159	5.8#	4.2cm	25m		56 box
JC0401510	6#	5.8cm	25m		32 box
JC0401511	7#	7.2cm	25m		32 box
JC0401512	8#	8.5cm	25m		32 box
JC0401513	9#	9.0cm	25m		24 box
JC0401514	10#	10cm	25m		24 box